AER BEATHA LTD

INNOVATIVE ASTHMA TREATMENT

Aer8™

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Effect of a Swedish-made user-friendly cardboard spacer on salbutamol MDI performance

Michael Friel MPharm PhD, Aer Beatha AB, Goteborg, Sweden

michael.friel@aerbeatha.com

Figure

Background

Spacers were developed to address problems of coordination and extrapulmonary drug deposition for patients using metered dose inhalers (MDIs). Respiratory experts recommend that every patient needing an MDI should regularly use a spacer, and know how to use it effectively, both for routine therapy and in an emergency (1).

Current plastic spacers are bulky, not very portable and need to be cleaned. Patients tend not to use them. The Aer8 spacer (Figure) addresses these shortcomings. Made in Sweden to Good Manufacturing Practice standards, it can be popped-up to erect and collapsed between uses. It is easily portable, large volume (500mL), non-static, compostable and much more environmentally friendly than plastic spacers.

Methods

Particle size distributions from six salbutamol (Ventolin®) MDIs were measured using a Next Generation Cascade Impactor® (NGI) set to an inlet flow rate of 30L/min and analyzed using high-performance liquid chromatography.

Based on the average mass of drug deposited in the actuator, spacer, mouth & throat and eight NGI stages, the Mass Median Aerodynamic Distribution (MMAD), Geometric Standard Deviation (GSD) and Fine Particle Dose (FPD, Respirable Dose) of the drug particles were calculated along with the Relative Standard Distribution percentage (%RSD).

Objective

•To compare drug particle size distribution and mouth and throat deposition from a salbutamol MDI with and without use of an Aer8 spacer.

Table of Results

Drug deposition area	Ventolin MDI alone (n=6)		Ventolin MDI + Aer8 Spacer (n=6)	
	Mean drug deposition (µg)	Mean %RSD	Mean drug deposition (µg)	Mean %RSD
Actuator	16	10	15	8
Spacer	Not Applicable		10	32
Mouth & Throat	48	22	2	30
NGI Stage 1	1	24	0	17
NGI Stage 2	1	25	1	13
NGI Stage 3	2	19	3	13
NGI Stage 4	12	20	17	12
NGI Stage 5	14	9	17	15
NGI Stage 6	4	8	4	12
NGI Stage 7	1	28	1	10
NGI Stage 8	1	26	0	5
Total deposition	100	5	69	11
MMAD µm	2	4	2	3
GSD	2	1	2	2
FPD µg (≤5.0µm)	32	13	40	11

Mean drug deposition (µg), MMAD (µm) GSD (µm) and FPD (µg) from salbutamol (Ventolin) MDIs alone and using an Aer8 spacer.

Interpretation of Results

- Compared with the salbutamol MDI used alone, use with an Aer8 spacer reduced mouth & throat (oropharyngeal) impaction, (mean deposition 48 µg v 2 µg), (Table).
- Total drug per actuation was reduced by the Aer8 Spacer, 69µg v 100µg with the MDI alone owing to the 2-second delay in sampling and aerosol depletion occurring with spacers (2).
- The MMAD, GSD and salbutamol respirable dose (FPD, mean 32µg v 40µg) were stable.

Conclusions

- Use of the Aer8 spacer with a salbutamol MDI reduced oropharyngeal particle impaction, with respirable dose delivery equivalent to the MDI alone.
- The Aer8 spacer respirable fine particle dose is similar to other spacers used to deliver salbutamol (3).
- With its user-friendly, environmentally-sustainable features, the Aer8 may be a viable alternative to existing plastic spacers and a useful tool for teaching patients how to use MDI devices effectively.

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